



5.0 검증 패널 3: 백혈구 활성화 (Leukocyte / Buffy Coat Activation)

사용자의 요구(모든 버피코트 활성화)는 L-PRP(백혈구-풍부 PRF) 또는 BMAC에 포함된 백혈구(Leukocytes: 단핵구, 림프구, 호중구)가 물리적 스트레스에 의해 활성화됨을 증명해야 함을 의미한다.




Scanning electron micrograph of a white blood cell (WBC) and platelets



Immunohistochemistry showing CD3 and CD20 expression within the fibrin matrix


A: CD 3 T cell



B: CD 20 B Cell

5.1 이론적 근거: 백혈구 기계 민감성 (Leukocyte Mechanosensitivity)

혈소판과 마찬가지로, 백혈구 역시 물리적 힘에 매우 민감한 세포이다.



2 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 4.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC


Noise Reduction = Line Int. Busy Gun Vacuum = 1.26e-009 mbar

Scan Speed = 3 N = 30

System Vacuum = 7.59e-007 mbar



Column Mode = Analytic



2 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 4.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.28e-009 mbar

Scan Speed = 3

N = 30

System Vacuum = 7.61e-007 mbar





2 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 4.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy Gun Vacuum = 1.27e-009 mbar

Scan Speed = 3 N = 30

System Vacuum = 7.62e-007 mbar



Column Mode = Analytic



2 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 5.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.27e-009 mbar


Scan Speed = 3

N = 30

System Vacuum = 7.62e-007 mbar



Column Mode = Analytic



1 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 5.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.26e-009 mbar


Scan Speed = 3

N = 30

System Vacuum = 7.56e-007 mbar



Column Mode = Analytic



2 μ m



EHT = 7.00 kV
WD = 14.2 mm
I Probe = 146 pA

Mag = 4.00 K X

Noise Reduction = Line Int. Busy

Scan Speed = 3

Signal A = SE2

N = 30

Date : 8 Jul 2025


Gun Vacuum = 1.18e-009 mbar

System Vacuum = 7.42e-007 mbar

YUMC



Column Mode = Analytic



2 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 4.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.28e-009 mbar

Scan Speed = 3

N = 30

System Vacuum = 7.61e-007 mbar



Column Mode = Analytic



2 μ m



EHT = 7.00 kV
WD = 11.4 mm
I Probe = 146 pA

Mag = 5.00 K X

Signal A = SE2

Date :8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.27e-009 mbar


ZEISS

Scan Speed = 3

N = 30

System Vacuum = 7.62e-007 mbar

Column Mode = Analytic



2 μ m



EHT = 7.00 kV
WD = 14.2 mm
I Probe = 146 pA

Mag = 4.00 K X

Signal A = SE2

Date : 8 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.25e-009 mbar


ZEISS

Scan Speed = 3

N = 30

System Vacuum = 7.53e-007 mbar

Column Mode = Analytic



1 μm



Column Mode = Analytic

EHT = 7.00 kV

WD = 14.2 mm

I Probe = 146 pA

Mag = 4.00 K X

Noise Reduction = Line Int. Busy

Scan Speed = 3

Signal A = SE2

N = 30


Date : 8 Jul 2025

Gun Vacuum = 1.22e-009 mbar

System Vacuum = 7.46e-007 mbar

YUMC





1 μ m

EHT = 5.00 kV

WD = 10.5 mm

I Probe = 146 pA

Mag = 7.00 K X

Noise Reduction = Line Int. Busy

Scan Speed = 3

Signal A = SE2

N = 30

Date :17 Jul 2025


Gun Vacuum = 1.12e-009 mbar

System Vacuum = 7.01e-007 mbar


YUMC

ZEISS

Column Mode = Analytic



1 μ m



EHT = 5.00 kV
WD = 10.5 mm
I Probe = 146 pA

Mag = 6.00 K X

Signal A = SE2

Date :17 Jul 2025

YUMC

Noise Reduction = Line Int. Busy

Gun Vacuum = 1.13e-009 mbar

Scan Speed = 3

N = 30

System Vacuum = 7.00e-007 mbar

Column Mode = Analytic

